



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



NEUROSURGERY

Neurosurgery deals with the diagnosis and treatment of pathological processes, which affect the nervous system. It includes the operative, non-operative, intensive care management and rehabilitation of patients with disorders affecting the brain and skull, spine and nervous system.

Perfection is in high demand in this discipline and requires talent, expertise and precision instrument manufacturing. we support the skilled hands of the surgeon with high quality instruments for perfect surgical work.

Our neurosurgery instruments include but are not limited to the following Special Micro- Neurosurgical Instruments, Special Microsurgical Tube- shaft Instruments, Head-Cranio Instruments and Holding Systems, Pituitary Instruments, Self- retaining Retractors, Spine Instruments, Rongeurs and Laminectomy Punches, Vessel- and Bulldog- Clamps, Vessel Clips, General Surgical Instruments, Neuro-Orthopedic Instruments and Headrest systems to fix the patient's head and flexible Leyla retractor systems.

Neurochirurgie

Neurocirugía

Neurochirurgie

Neurochirurgia



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Self-Retaining Brain Retractor
Selbsthaltender Hirnretractor
Separador cerebral autostático
Retracteur cérébrale autostatique
Retratore celebrale autostatico

Yasargil (Leyla) Brain Retractor Single Arm Complete

Catalog #	Size
NS-010-01	N/A




Yasargil (Leyla) Brain Double Retractor Arm Complete

Catalog #	Size
NS-010-02	N/A



Yasargil (Leyla) Flexible Arm

Catalog #	Size
NS-010-10	N/A



Yasargil (Leyla) Flexible Base


Catalog #	Size
NS-010-21	N/A
NS-010-22	N/A
NS-010-30	N/A



NS-010-21



NS-010-22



NS-010-30

Yasargil (Leyla) Support for Round Spatulas

Catalog #	Size
NS-010-40	N/A





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Yasargil (Leyla) Support for Flat Spatulas

Catalog #	Size
NS-010-45	N/A

Self-Retaining Brain Retractor
Selbsthaltender Hirnretractor
Separador cerebral autostatico
Retracteur cérébrale autostatique
Retrattore celebrale autostatico



Yasargil (Leyla) Ball and Socket Joint

Catalog #	Size
NS-012-09	9 x 32 mm
NS-012-10	10 x 25 mm



Yasargil (Leyla) Holding Rod

Catalog #	Size
NS-012-20	N/A



Yasargil (Leyla) Coupling Head

Catalog #	Size
NS-012-25	N/A
NS-012-40	N/A



NS-012-25



NS-012-40

Caspar Vertebral Distractor

Right

Catalog #	Size
NS-020-01	N/A



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Vertebral Distractor
Wirbel distractor
Vertebrales distractor
Vertébrale distracteur
Vertebrale distrattore

Caspar Drill Guide

Right

Catalog #	Size
NS-020-02	N/A



Casper Vertebral Distractor

Left

Catalog #	Size
NS-022-01	N/A



Caspar Drill Guide

Left

Catalog #	Size
NS-022-02	N/A




Caspar Twist Drill

Catalog #	Size
NS-024-00	14.5 cm (5 1/2")



Caspar Distraction Screw

Catalog #	Working End
NS-026-14	14 mm
NS-026-16	16 mm
NS-026-18	18 mm





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Casper Screw Driver


Wire Saws, Skull Traction Tongs
Drahtsägen, Schädelhalter
Sierras de hilo, Estribos de tracción craneal
Fils-scies, Pinces de traction crânienne
Fili sega, Pinze per trazione craniale



Catalog #	Size
NS-028-20	20 cm (8")

Gigli Wire Saw

Catalog #	Size
NS-050-30	30 cm (11 3/4")
NS-050-40	40 cm (11 3/4")
NS-050-50	50 cm (11 3/4")
NS-050-60	60 cm (11 3/4")
NS-050-70	70 cm (11 3/4")



Gigli Wire Saw Handle

Catalog #	Size
NS-055-00	N/A



De Martel Wire Saw Guide



Catalog #	Size
NS-058-00	33 cm (13")

Crutchfield Skull Traction Tongs

Small Pattern

Catalog #	Size
NS-061-01	N/A



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.

Manufacturers of single use and reusable instruments



Skull Traction Tongs

Schädelhalter

Estribos de tracción craneal

Pinces de traction crânienne

Pinze per trazione craniale

Crutchfield Skull Traction Tongs

Medium Pattern

Catalog #	Size
NS-061-02	N/A



Crutchfield Drill Point

Catalog #	Size, Diameter
NS-065-35	7 cm (2 3/4"), Ø 3.5 mm
NS-065-50	7 cm (2 3/4"), Ø 5.0 mm



Skull Traction Tongs

Large Pattern

Catalog #	Size
NS-071-00	N/A



Spare Pin

Catalog #	Size
NS-075-01	N/A



Wrench

Catalog #	Size
NS-075-02	N/A





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Dahlgren Cranial Rongeur

Cranial Rongeurs
Schädelzangen
Pinzas sacabocados craneales
Pincas emporte-pièces crâniennes
Pinze per craniectomia



Catalog #	Size
NS-091-20	20 cm (8")

De Vilbiss Cranial Rongeur



Catalog #	Size
NS-093-21	20.5 cm (8")

Cone Cranial Punch



Catalog #	Size
NS-095-22	22.5 cm (8 3/4")

Raney Scalp Clip



Catalog #	Size
NS-100-00	N/A

Raney Applying - Removal Forcep



Catalog #	Size
NS-102-00	15.5 cm (6")

NEUROSURGERY




CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Clips and Applying Forceps
Clips und Anlegezangen
Grapas y pinzas aplicadoras
Clips et pinces pour appliquer
Agrafes e pinze per applicare

McKenzie Clip Applying Forcep



Catalog #	Size
NS-112-19	18.5 cm (7 1/4")

Olivecrona-Toennis Clip Applying Forcep



Catalog #	Size, Shaft Length
NS-115-14	N/A, 140 mm


Silver Brain Clip

Pack of 100




Catalog #	Size
NS-118-00	N/A

Yasargil Galea Hook



Catalog #	Size
NS-120-01	N/A
NS-120-02	N/A

Davis Dura Dissector



Catalog #	Size
NS-151-21	21.5 cm (8 1/2")



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Davis Dura Dissector

Dura Dissectors
Duradisektoren
Disecadores para duramadre
Dissecteurs à dure-mère
Dissettori per la dura-madre




Catalog #	Size
NS-153-24	24.5 cm (9 ¾")

Frazier Dura Dissector




Catalog #	Size
NS-155-01	15 cm (6")

Hoehn Dura Dissector



NS-155-02



NS-155-03



Catalog #	Size
NS-155-02	15 cm (6")
NS-155-03	15 cm (6")

Love-Adson Dura Dissector



Catalog #	Size
NS-157-16	16 cm (6 ½")

Toennis Dura Dissector



Catalog #	Size
NS-159-26	26 cm (10 ¼")

NEUROSURGERY




CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Dura Dissectors
Duradisektoren
Disecadores para duramadre
Dissecteurs à dure-mère
Dissettori per la dura-madre

Olivecrona Dura Dissector



Catalog #	Size
NS-161-01	24 cm (9 1/2")

Olivecrona Dura Dissector



Catalog #	Size
NS-161-02	18 cm (7")

Olivecrona Dura Dissector



Catalog #	Size
NS-161-03	24 cm (9 1/2")

Olivecrona Dura Dissector



Catalog #	Size
NS-161-04	20 cm (8")

Sachs Dura Dissector

Sharp



Catalog #	Size
NS-163-21	20 cm (8")



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Milligan Dura Dissector

Dura Dissectors
Duradissektoren
Disecadores para duramadre
Dissecteurs à dure-mère
Dissettori per la dura-madre



Catalog #	Size
NS-165-21	21.5 cm (8 1/2")

Dura Dissector



Catalog #	Size
NS-167-22	22 cm (8 3/4")

Penfield Dura Dissector



Catalog #	Size, Profile
NS-171-01	17.5 cm (7"), Fig. 1

Penfield Dura Dissector



Catalog #	Size, Profile
NS-171-02	19.5 cm (7 3/4"), Fig. 2

Penfield Dura Dissector



Catalog #	Size, Profile
NS-171-03	19.5 cm (7 3/4"), Fig. 3

NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Dura Dissectors
Duradissektoren
Disecadores para duramadre
Dissecteurs à dure-mère
Dissettori per la dura-madre

Penfield Dura Dissector



Catalog #	Size, Profile
NS-171-04	20.5 cm (8"), Fig. 4

Penfield Dura Dissector



Catalog #	Size, Profile
NS-171-05	18 cm (7"), Fig. 5

Horsley Dura Elevator



Catalog #	Size
NS-181-20	20 cm (8")

Sachs Dura Elevator



Catalog #	Size
NS-183-21	21 cm (8 1/4")

Love Nerve Root Retractor

Straight



Catalog #	Size
NS-190-21	21.5 cm (8 1/2")



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Nerve Root Retractors
Nervwurzelhaken
Separadores para la raíz de nervio
Retracteurs pour racines nerveuses
Uncini per le radici die nervi

Love Nerve Root Retractor

Angled 45°


Catalog #	Size
NS-191-18	17.5 cm (7")



Love Nerve Root Retractor

Angled 90°

Catalog #	Size
NS-193-11	11 cm (4 1/4")



Love Nerve Root Retractor

Catalog #	Size
NS-195-21	22 cm (8 3/4")



Hess Nerve Root Retractor

Catalog #	Size, Width
NS-197-03	17.5 cm (7"), 3 mm
NS-197-05	17.5 cm (7"), 5 mm
NS-197-07	17.5 cm (7"), 7 mm



Scoville Nerve Root Retractor

Catalog #	Size
NS-200-23	23.5 cm (9 1/4")



NEUROSURGERY



Nerve Root Retractors
Nervwurzelnhaken
Separadores para la raíz de nervio
Retracteurs pour racines nerveuses
Uncini per le radici die nervi

Holsher Nerve Root Retractor



Catalog #	Size
NS-201-24	23.5 cm (9 1/4")

Scoville Nerve Root Retractor



Catalog #	Size
NS-203-19	16.5 cm (6 1/2")


Cushing Brain Spatulas



Catalog #	Size, Blade Size
NS-210-07	18 cm (7"), 7 + 9 mm
NS-210-11	18 cm (7"), 11 + 13 mm
NS-210-15	18 cm (7"), 15 + 18 mm
NS-210-18	18 cm (7"), 18 + 22 mm

Olivecrona Brain Spatulas

Concave



Catalog #	Size, Blade Size
NS-213-07	18 cm (7"), 7 + 9 mm
NS-213-11	18 cm (7"), 11 + 13 mm
NS-213-15	18 cm (7"), 15 + 18 mm
NS-213-18	18 cm (7"), 18 + 22 mm



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Olivecrona Brain Spatulas

Convex

Brain Spatulas
Hirnspatel
Espátulas cerebrales
Spatules cérébrales
Spatole cerebrali

Catalog #	Size, Blade Size
NS-215-07	18 cm (7"), 7 + 9 mm
NS-215-11	18 cm (7"), 11 + 13 mm
NS-215-15	18 cm (7"), 15 + 18 mm
NS-215-18	18 cm (7"), 18 + 22 mm



Heifetz Brain Spatulas

Catalog #	Size, Blade Size
NS-218-07	20 cm (8"), 7 + 8 mm
NS-218-10	20 cm (8"), 10 + 11 mm
NS-218-13	20 cm (8"), 13 + 14 mm
NS-218-15	20 cm (8"), 15 + 17 mm
NS-218-19	20 cm (8"), 19 + 20 mm
NS-218-22	20 cm (8"), 22 + 22 mm



Heifetz Brain Spatulas

Round Handle

Catalog #	Size, Blade Size
NS-221-08	20 cm (8"), 8 mm
NS-221-11	20 cm (8"), 11 mm
NS-221-14	20 cm (8"), 14 mm
NS-221-17	20 cm (8"), 17 mm
NS-221-20	20 cm (8"), 20 mm



Brain Spatulas

Malleable

Catalog #	Size
NS-224-12	20 cm (8"), 12 mm
NS-224-17	20 cm (8"), 17 mm



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Nerve Hooks
Nervhaken
Ganchos para nervios
Crochets à nerfs
Uncini per nervi

Smithwick Nerve Hook




Catalog #	Size
NS-231-17	17 cm (6 ¾")

Smithwick Nerve Hook



Catalog #	Size
NS-233-21	21 cm (8 ¼")

Adson Dura Hook



Catalog #	Size
NS-235-20	20 cm (8")

Graham Dura Hook

Blunt



Catalog #	Size
NS-237-16	16.5 cm (6 ½")

Sachs Dura Hook

Sharp



Catalog #	Size
NS-239-18	18 cm (7")



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Frazier Dura Hook

Sharp

Nerve Hooks
Nervhaken
Ganchos para nervios
Crochets à nerfs
Uncini per nervi



Catalog #	Size
NS-241-18	18 cm (7")

Hoehn Nerve Retractor



Catalog #	Size
NS-245-22	22 cm (8 3/4")

Dandy Nerve Retractor



NS-251-01



NS-251-02



NS-251-03

Catalog #	Size
NS-251-01	20 cm (8")
NS-251-02	20 cm (8")
NS-251-03	20 cm (8")

Smithwick Nerve Retractor



Catalog #	Size
NS-253-29	29 cm (11 1/2")

Toennis Nerve Retractor



Catalog #	Size
NS-261-15	15 cm (6")

NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Dura Knives, Laminectomy Punches
Duramesser, Laminektomiestanzen
Bisturís para duramadre, Pinzas cortantes para laminectomía
Bistouries pour duremère, Pinces emporte-pièces à laminectomie
Bisturi per duramadre, Pinze taglianti per laminectomia

Toennis Brain Knife



Catalog #	Size
NS-280-18	18 cm (7")


Halle Brain Knife



Catalog #	Size
NS-282-20	20 cm (8")

Ferris-Smith Kerrison Rongeur


Up Cutting



Catalog #	Size, Bite Size
NS-300-01	15 cm (6"), 1 mm
NS-300-02	15 cm (6"), 2 mm
NS-300-03	15 cm (6"), 3 mm
NS-300-04	15 cm (6"), 4 mm
NS-300-05	15 cm (6"), 5 mm

Ferris-Smith Kerrison Rongeur

Down Cutting



Catalog #	Size, Bite Size
NS-301-01	15 cm (6"), 1 mm
NS-301-02	15 cm (6"), 2 mm
NS-301-03	15 cm (6"), 3 mm
NS-301-04	15 cm (6"), 4 mm
NS-301-05	15 cm (6"), 5 mm




CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Ferris-Smith Kerrison Rongeur

40° Forward Up Cutting


Catalog #	Size, Bite Size
NS-302-01	15 cm (6"), 1 mm
NS-302-02	15 cm (6"), 2 mm
NS-302-03	15 cm (6"), 3 mm
NS-302-04	15 cm (6"), 4 mm
NS-302-05	15 cm (6"), 5 mm



Ferris-Smith Kerrison Rongeur

40° Forward Down Cutting


Catalog #	Size, Bite Size
NS-303-01	15 cm (6"), 1 mm
NS-303-02	15 cm (6"), 2 mm
NS-303-03	15 cm (6"), 3 mm
NS-303-04	15 cm (6"), 4 mm
NS-303-05	15 cm (6"), 5 mm



Ferris-Smith Kerrison Rongeur

Up Cutting

Catalog #	Size, Bite Size
NS-304-01	18 cm (7"), 1 mm
NS-304-02	18 cm (7"), 2 mm
NS-304-03	18 cm (7"), 3 mm
NS-304-04	18 cm (7"), 4 mm
NS-304-05	18 cm (7"), 5 mm



NEUROSURGERY




CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pincès emporte-pièces à laminectomie
Pinze taglianti per laminectomia


Ferris-Smith Kerrison Rongeur
Down Cutting

Catalog #	Size, Bite Size
NS-305-01	18 cm (7"), 1 mm
NS-305-02	18 cm (7"), 2 mm
NS-305-03	18 cm (7"), 3 mm
NS-305-04	18 cm (7"), 4 mm
NS-305-05	18 cm (7"), 5 mm




Ferris-Smith Kerrison Rongeur
40° Forward Up Cutting

Catalog #	Size, Bite Size
NS-306-01	18 cm (7"), 1 mm
NS-306-02	18 cm (7"), 2 mm
NS-306-03	18 cm (7"), 3 mm
NS-306-04	18 cm (7"), 4 mm
NS-306-05	18 cm (7"), 5 mm



Ferris-Smith Kerrison Rongeur
40° Forward Down Cutting

Catalog #	Size, Bite Size
NS-307-01	18 cm (7"), 1 mm
NS-307-02	18 cm (7"), 2 mm
NS-307-03	18 cm (7"), 3 mm
NS-307-04	18 cm (7"), 4 mm
NS-307-05	18 cm (7"), 5 mm






CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Ferris-Smith Kerrison Rongeur

Up Cutting

Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pinces emporte-pièces à laminectomie
Pinze taglianti per laminectomia




Catalog #	Size, Bite Size
NS-308-01	20 cm (8"), 1 mm
NS-308-02	20 cm (8"), 2 mm
NS-308-03	20 cm (8"), 3 mm
NS-308-04	20 cm (8"), 4 mm
NS-308-05	20 cm (8"), 5 mm




Ferris-Smith Kerrison Rongeur

Down Cutting




Catalog #	Size, Bite Size
NS-309-01	20 cm (8"), 1 mm
NS-309-02	20 cm (8"), 2 mm
NS-309-03	20 cm (8"), 3 mm
NS-309-04	20 cm (8"), 4 mm
NS-309-05	20 cm (8"), 5 mm




Ferris-Smith Kerrison Rongeur

40° Forward Up Cutting



Catalog #	Size, Bite Size
NS-310-01	20 cm (8"), 1 mm
NS-310-02	20 cm (8"), 2 mm
NS-310-03	20 cm (8"), 3 mm
NS-310-04	20 cm (8"), 4 mm
NS-310-05	20 cm (8"), 5 mm



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.

Manufacturers of single use and reusable instruments




Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pincers emporte-pièces à laminectomie
Pinze taglienti per laminectomia

Ferris-Smith Kerrison Rongeur

40° Forward Down Cutting

Catalog #	Size, Bite Size
NS-311-01	20 cm (8"), 1 mm
NS-311-02	20 cm (8"), 2 mm
NS-311-03	20 cm (8"), 3 mm
NS-311-04	20 cm (8"), 4 mm
NS-311-05	20 cm (8"), 5 mm



Ferris-Smith Kerrison Rongeur

Detachable Model, Up Cutting

Catalog #	Size, Bite Size
NS-320-01	18 cm (7"), 1 mm
NS-320-02	18 cm (7"), 2 mm
NS-320-03	18 cm (7"), 3 mm
NS-320-04	18 cm (7"), 4 mm
NS-320-05	18 cm (7"), 5 mm



Ferris-Smith Kerrison Rongeur

Detachable Model, Down Cutting

Catalog #	Size, Bite Size
NS-321-01	18 cm (7"), 1 mm
NS-321-02	18 cm (7"), 2 mm
NS-321-03	18 cm (7"), 3 mm
NS-321-04	18 cm (7"), 4 mm
NS-321-05	18 cm (7"), 5 mm





CENTRAL SURGICAL CO.(PVT)LTD.

Manufacturers of single use and reusable instruments



Laminectomy Punches

Laminektomiestanzen

Pinzas cortantes para laminectomía

Pinces emporte-pièces à laminectomie

Pinze taglianti per laminectomia

Ferris-Smith Kerrison Rongeur

Detachable Model, 40° Forward Up Cutting

Catalog

Size, Bite Size

NS-322-01	18 cm (7"), 1 mm
NS-322-02	18 cm (7"), 2 mm
NS-322-03	18 cm (7"), 3 mm
NS-322-04	18 cm (7"), 4 mm
NS-322-05	18 cm (7"), 5 mm



Ferris-Smith Kerrison Rongeur

Detachable Model, 40° Forward Down Cutting

Catalog

Size, Bite Size

NS-323-01	18 cm (7"), 1 mm
NS-323-02	18 cm (7"), 2 mm
NS-323-03	18 cm (7"), 3 mm
NS-323-04	18 cm (7"), 4 mm
NS-323-05	18 cm (7"), 5 mm



Ferris-Smith Kerrison Rongeur

Detachable Model, Up Cutting

Catalog

Size, Bite Size

NS-324-01	20 cm (8"), 1 mm
NS-324-02	20 cm (8"), 2 mm
NS-324-03	20 cm (8"), 3 mm
NS-324-04	20 cm (8"), 4 mm
NS-324-05	20 cm (8"), 5 mm



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.

Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pincès emporte-pièces à laminectomie
Pinze taglienti per laminectomia

Ferris-Smith Kerrison Rongeur

Detachable Model, Down Cutting

Catalog #	Size, Bite Size
NS-325-01	20 cm (8"), 1 mm
NS-325-02	20 cm (8"), 2 mm
NS-325-03	20 cm (8"), 3 mm
NS-325-04	20 cm (8"), 4 mm
NS-325-05	20 cm (8"), 5 mm



Ferris-Smith Kerrison Rongeur

Detachable Model, 40° Up Cutting

Catalog #	Size, Bite Size
NS-326-01	20 cm (8"), 1 mm
NS-326-02	20 cm (8"), 2 mm
NS-326-03	20 cm (8"), 3 mm
NS-326-04	20 cm (8"), 4 mm
NS-326-05	20 cm (8"), 5 mm



Ferris-Smith Kerrison Rongeur

Detachable Model, 40° Down Cutting

Catalog #	Size, Bite Size
NS-327-01	20 cm (8"), 1 mm
NS-327-02	20 cm (8"), 2 mm
NS-327-03	20 cm (8"), 3 mm
NS-327-04	20 cm (8"), 4 mm
NS-327-05	20 cm (8"), 5 mm





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pincas emporte-pièces à laminectomie
Pinze taglienti per laminectomia

1 mm 2 mm 3 mm 4 mm 5 mm

Ferris-Smith Kerrison Rongeur

Up Cutting



Catalog

Size, Bite Size

NS-330-01	18 cm (7"), 1 mm
NS-330-02	18 cm (7"), 2 mm
NS-330-03	18 cm (7"), 3 mm
NS-330-04	18 cm (7"), 4 mm
NS-330-05	18 cm (7"), 5 mm

Ferris-Smith Kerrison Rongeur

Down Cutting

1 mm 2 mm 3 mm 4 mm 5 mm



Catalog

Size, Bite Size

NS-331-01	18 cm (7"), 1 mm
NS-331-02	18 cm (7"), 2 mm
NS-331-03	18 cm (7"), 3 mm
NS-331-04	18 cm (7"), 4 mm
NS-331-05	18 cm (7"), 5 mm

Ferris-Smith Kerrison Rongeur

40° Forward Up Cutting

1 mm 2 mm 3 mm 4 mm 5 mm



Catalog

Size, Bite Size

NS-332-01	18 cm (7"), 1 mm
NS-332-02	18 cm (7"), 2 mm
NS-332-03	18 cm (7"), 3 mm
NS-332-04	18 cm (7"), 4 mm
NS-332-05	18 cm (7"), 5 mm

NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pinces emporte-pièces à laminectomie
Pinze taglianti per laminectomia

Ferris-Smith Kerrison Rongeur

40° Forward Down Cutting


Catalog #	Size, Bite Size
NS-333-01	18 cm (7"), 1 mm
NS-333-02	18 cm (7"), 2 mm
NS-333-03	18 cm (7"), 3 mm
NS-333-04	18 cm (7"), 4 mm
NS-333-05	18 cm (7"), 5 mm



Ferris-Smith Kerrison Rongeur

Up Cutting


Catalog #	Size, Bite Size
NS-334-01	20 cm (8"), 1 mm
NS-334-02	20 cm (8"), 2 mm
NS-334-03	20 cm (8"), 3 mm
NS-334-04	20 cm (8"), 4 mm
NS-334-05	20 cm (8"), 5 mm



Ferris-Smith Kerrison Rongeur

Down Cutting

Catalog #	Size, Bite Size
NS-335-01	20 cm (8"), 1 mm
NS-335-02	20 cm (8"), 2 mm
NS-335-03	20 cm (8"), 3 mm
NS-335-04	20 cm (8"), 4 mm
NS-335-05	20 cm (8"), 5 mm





CENTRAL SURGICAL CO.(PVT)LTD.

Manufacturers of single use and reusable instruments



Ferris-Smith Kerrison Rongeur

40° Forward Up Cutting


Laminectomy Punches

Laminektomiestanzen

Pinzas cortantes para laminectomía

Pinces emporte-pièces à laminectomie

Pinze taglienti per laminectomia




Catalog

Size, Bite Size

NS-336-01	20 cm (8"), 1 mm
NS-336-02	20 cm (8"), 2 mm
NS-336-03	20 cm (8"), 3 mm
NS-336-04	20 cm (8"), 4 mm
NS-336-05	20 cm (8"), 5 mm

Ferris-Smith Kerrison Rongeur

40° Forward Down Cutting



Catalog

Size, Bite Size

NS-337-01	20 cm (8"), 1 mm
NS-337-02	20 cm (8"), 2 mm
NS-337-03	20 cm (8"), 3 mm
NS-337-04	20 cm (8"), 4 mm
NS-337-05	20 cm (8"), 5 mm

Cushing Leminectomy Rongeur

Catalog

Size, Bite Size, Profile

NS-400-01	13 cm (5"), 2 x 10 mm, Straight
NS-400-02	13 cm (5"), 2 x 10 mm, Up
NS-400-03	13 cm (5"), 2 x 10 mm, Down



Cushing Leminectomy Rongeur

Catalog

Size, Bite Size, Profile

NS-401-01	15 cm (6"), 2 x 10 mm, Straight
NS-401-02	15 cm (6"), 2 x 10 mm, Up
NS-401-03	15 cm (6"), 2 x 10 mm, Down



NEUROSURGERY




CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pinces emporte-pièces à laminectomie
Pinze taglienti per laminectomia


Cushing Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-402-01	18 cm (7"), 2 x 10 mm, Straight
NS-402-02	18 cm (7"), 2 x 10 mm, Up
NS-402-03	18 cm (7"), 2 x 10 mm, Down



Cushing Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-403-01	20 cm (8"), 2 x 10 mm, Straight
NS-403-02	20 cm (8"), 2 x 10 mm, Up
NS-403-03	20 cm (8"), 2 x 10 mm, Down





Love-Gruenwald Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-404-01	13 cm (5"), 3 x 10 mm, Straight
NS-404-02	13 cm (5"), 3 x 10 mm, Up
NS-404-03	13 cm (5"), 3 x 10 mm, Down



Love-Gruenwald Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-405-01	15 cm (6"), 3 x 10 mm, Straight
NS-405-02	15 cm (6"), 3 x 10 mm, up
NS-405-03	15 cm (6"), 3 x 10 mm, Down





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pinces emporte-pièces à laminectomie
Pinze taglienti per laminectomia

Love-Gruenwald Laminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-406-01	18 cm (7"), 3 x 10 mm, Straight
NS-406-02	18 cm (7"), 3 x 10 mm, Up
NS-406-03	18 cm (7"), 3 x 10 mm, Down



Love-Gruenwald Laminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-407-01	20 cm (8"), 3 x 10 mm, Straight
NS-407-02	20 cm (8"), 3 x 10 mm, Up
NS-407-03	20 cm (8"), 3 x 10 mm, Down




Spurling Laminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-408-01	13 cm (5"), 4 x 10 mm, Straight
NS-408-02	13 cm (5"), 4 x 10 mm, Up
NS-408-03	13 cm (5"), 4 x 10 mm, Down



Spurling Laminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-409-01	15 cm (6"), 4 x 10 mm, Straight
NS-409-02	15 cm (6"), 4 x 10 mm, Up
NS-409-03	15 cm (6"), 4 x 10 mm, Down



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pinces emporte-pièces à laminectomie
Pinze taglienti per laminectomia

Spurling Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-410-01	18 cm (7"), 4 x 10 mm, Straight
NS-410-02	18 cm (7"), 4 x 10 mm, up
NS-410-03	18 cm (7"), 4 x 10 mm, Down




Spurling Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-411-01	20 cm (8"), 4 x 10 mm, Straight
NS-411-02	20 cm (8"), 4 x 10 mm, Up
NS-411-03	20 cm (8"), 4 x 10 mm, Down




Ferris-Smith Leminectomy Rongeur Straight

Catalog #	Size, Bite Size, Profile
NS-430-01	15.5 cm (6"), 2 mm, Straight
NS-430-02	15.5 cm (6"), 2 mm, Up
NS-430-03	15.5 cm (6"), 2 mm, Down



Ferris-Smith Leminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-431-01	15.5 cm (6"), 3 mm, Straight
NS-431-02	15.5 cm (6"), 3 mm, Up
NS-431-03	15.5 cm (6"), 3 mm, Down





CENTRAL SURGICAL CO.(PVT)LTD.


Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pinces emporte-pièces à laminectomie
Pinze taglienti per laminectomia


Ferris-Smith Laminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-432-01	15.5 cm (6"), 4 mm, Straight
NS-432-02	15.5 cm (6"), 4 mm, Up
NS-432-03	15.5 cm (6"), 4 mm, Down



Ferris-Smith Laminectomy Rongeur

Catalog #	Size, Bite Size, Profile
NS-433-01	15.5 cm (6"), 5 mm, Straight
NS-433-02	15.5 cm (6"), 5 mm, Up
NS-433-03	15.5 cm (6"), 5 mm, Down



Laminectomy Rongeur

Fenestrated and Serrated Jaws

Catalog #	Size, Bite Size, Profile
NS-440-01	15.5 cm (6"), 2 x 12 mm, Straight
NS-440-02	15.5 cm (6"), 2 x 12 mm, Up
NS-440-03	15.5 cm (6"), 2 x 12 mm, Down



Laminectomy Rongeur

Fenestrated and Serrated Jaws

Catalog #	Size, Bite Size, Profile
NS-441-01	15.5 cm (6"), 3 x 12 mm, Straight
NS-441-02	15.5 cm (6"), 3 x 12 mm, Up
NS-441-03	15.5 cm (6"), 3 x 12 mm, Down



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Laminectomy Punches
Laminektomiestanzen
Pinzas cortantes para laminectomía
Pincès emporte-pièces à laminectomie
Pinze taglienti per laminectomia

Laminectomy Rongeur

Straight, Fenestrated and Serrated Jaws

Catalog #	Size, Bite Size, Profile
NS-442-01	15.5 cm (6"), 4 x 14 mm, Straight
NS-442-02	15.5 cm (6"), 4 x 14 mm, Up
NS-442-03	15.5 cm (6"), 4 x 14 mm, Down



Laminectomy Rongeur

Fenestrated and Serrated Jaws

Catalog #	Size, Bite Size, Profile
NS-443-01	15.5 cm (6"), 5 x 14 mm, Straight
NS-443-02	15.5 cm (6"), 5 x 14 mm, Up
NS-443-03	15.5 cm (6"), 5 x 14 mm, Down



Laminectomy Rongeur

Fenestrated and Serrated Jaws

Catalog #	Size, Bite Size, Profile
NS-444-01	15.5 cm (6"), 6 x 14 mm, Straight
NS-444-02	15.5 cm (6"), 6 x 14 mm, Up
NS-444-03	15.5 cm (6"), 6 x 14 mm, Down



Schlesinger Laminectomy Rongeur

Serrated Jaws

Catalog #	Size, Bite Size
NS-451-02	13 cm (5"), 2 x 10 mm
NS-451-03	13 cm (5"), 3 x 10 mm
NS-451-04	13 cm (5"), 4 x 10 mm





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Yasargil Tumor Grasping Forcep


Laminectomy Punches, Micro Neuro Instruments
Laminektomiestanzen, Mikro Neuro Instrumente
Pinzas cortantes para laminectomía, Micro Neuro Instrumentos
Pincas emporte-pièces à laminectomie, Micro Neuro Instruments
Pinze taglianti per laminectomia, Micro Neuro Strumenti

Catalog #	Size, Bite Size
NS-461-35	20 cm (8"), Ø 3.5 mm



Oldberg Tumor Grasping Forcep

Catalog #	Size, Bite Size
NS-463-60	18 cm (7"), Ø 6.5 mm



Jacobson Micro Knife


Straight



Catalog #	Size
NS-550-18	18.5 cm (7 1/4")

Koos Micro Knife

Bayonet Shaped




Catalog #	Size, Bite Size
NS-555-01	18.5 cm (7 1/4")
NS-555-02	18.5 cm (7 1/4")
NS-555-03	18.5 cm (7 1/4")
NS-555-04	18.5 cm (7 1/4")



Koos Micro Knife

Bayonet Shaped



Catalog #	Size
NS-557-01	18.5 cm (7 1/4")
NS-557-02	18.5 cm (7 1/4")



NEUROSURGERY





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti




Yasargil Micro Knife
Bayonet Shaped



Catalog #	Size
NS-559-01	18.5 cm (7 1/4")
NS-559-02	18.5 cm (7 1/4")




Yasargil Micro Respiratory
Straight



Catalog #	Size
NS-564-01	18.5 cm (7 1/4")
NS-564-02	18.5 cm (7 1/4")
NS-564-03	18.5 cm (7 1/4")
NS-564-04	18.5 cm (7 1/4")
NS-564-05	18.5 cm (7 1/4")



Yasargil Micro Respiratory
Bayonet Shaped



Catalog #	Size
NS-565-01	18.5 cm (7 1/4")
NS-565-02	18.5 cm (7 1/4")
NS-565-03	18.5 cm (7 1/4")



Yasargil Micro Curette
Bayonet Shaped



Catalog #	Size
NS-567-18	18.5 cm (7 1/4")



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Micro Nerve & Vessel Hook



Catalog #	Size, Tip Length
NS-571-01	18.5 cm (7 1/4"), 3.0 mm
NS-571-02	18.5 cm (7 1/4"), 5.0 mm


Krayenbuhl Micro Nerve & Vessel Hook



Catalog #	Size, Profile
NS-573-01	18.5 cm (7 1/4"), Fig. 1
NS-573-02	18.5 cm (7 1/4"), Fig. 2
NS-573-03	18.5 cm (7 1/4"), Fig. 3
NS-573-04	18.5 cm (7 1/4"), Fig. 4



Yasargil Ligature Guide



Catalog #	Size, Profile
NS-575-01	18.5 cm (7 1/4"), Fig. 1
NS-575-02	18.5 cm (7 1/4"), Fig. 2
NS-575-03	18.5 cm (7 1/4"), Fig. 3
NS-575-04	18.5 cm (7 1/4"), Fig. 4



Micro Vessel Dilator



Catalog #	Size, Diameter, Profile
NS-577-01	12.5 cm (5"), Ø 0.5 mm, Fig. 1
NS-577-02	12.5 cm (5"), Ø 1.0 mm, Fig. 2
NS-577-3	12.5 cm (5"), Ø 1.5 mm, Fig. 3

NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.

Manufacturers of single use and reusable instruments



Micro Neuro Instruments


Mikro Neuro Instrumente

Micro Neuro Instrumentos

Micro Neuro Instruments

Micro Neuro Strumenti

Jacobson Micro Vessel Probe




Catalog


Size, Profile

NS-579-01 18.5 cm (7 1/4"), 1 Prong

NS-579-02 18.5 cm (7 1/4"), 2 Prong




Jacobson Micro Suture Pusher



Catalog

Size

NS-579-03 18.5 cm (7 1/4")



Rhoton Micro Dissector

Set of 15 Ref: NS-591-01 up to NS-591-15




Catalog

Size


NS-590-00 N/A

Rhoton Micro Dissector


Round Shaped



NS-591-01



NS-591-02



NS-591-03

Catalog

Size, Diameter

NS-591-01 18.5 cm (7 1/4"), Ø 1.0 mm


NS-591-02 18.5 cm (7 1/4"), Ø 2.0 mm

NS-591-03 18.5 cm (7 1/4"), Ø 3.0 mm




Rhoton Micro Dissector

Spatula Shaped



NS-591-04



NS-591-05



NS-591-06

Catalog

Size, Tip Size

NS-591-04 18.5 cm (7 1/4"), 1.0 mm

NS-591-05 18.5 cm (7 1/4"), 1.5 mm


NS-591-06 18.5 cm (7 1/4"), 2.0 mm





Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Rhoton Micro Elevator



Catalog #	Size, Tip Size
NS-591-07	18.5 cm (7 1/4"), 1.2 mm
NS-591-08	18.5 cm (7 1/4"), 2.8 mm

Rhoton Micro Needle

Semi, Sharp



Catalog #	Size
NS-591-09	18.5 cm (7 1/4")

Rhoton Micro Hook

Angled 90°, Blunt



Catalog #	Size, Tip Size
NS-591-10	18.5 cm (7 1/4"), 2 mm

Rhoton Micro Hook

Angled 90°, Semi Sharp



Catalog #	Size, Tip Size
NS-591-11	18.5 cm (7 1/4"), 2 mm

Rhoton Micro Hook

Angled 45°, Semi Sharp



Catalog #	Size, Tip Size
NS-591-12	18.5 cm (7 1/4"), 3 mm



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Rhoton Micro Hook

Angled 90°




Catalog #

Size

NS-591-13

18.5 cm (7 1/4")

Rhoton Micro Curette



Catalog #


Size, Tip Size

NS-591-14

18.5 cm (7 1/4"), 1.0 x 2.0 mm

Rhoton Micro Curette

Angled 45°



Catalog #


Size, Tip Size

NS-591-15

18.5 cm (7 1/4"), 1.0 x 2.0 mm

Yasargil Micro Elevator

Bayonet Shaped



Catalog #


Size

NS-607-24

24 cm (9 1/2")

Yasargil Micro Elevator

Bayonet Shaped



Catalog #

Size

NS-611-24

24 cm (9 1/2")



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Yasargil Micro Raspatory

Bayonet Shaped

Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Catalog #	Size
NS-613-24	24 cm (9 1/2")



Yasargil Micro Dissector

Bayonet Shaped

Catalog #	Size
NS-615-01	24 cm (9 1/2")
NS-615-02	24 cm (9 1/2")



Yasargil Micro Needle

Bayonet Shaped, Semi Sharp

Catalog #	Size
NS-618-24	24 cm (9 1/2")



Yasargil Micro Curette

Straight

Catalog #	Size
NS-631-01	24 cm (9 1/2")



Yasargil Micro Curette

Angled

Catalog #	Size
NS-633-01	24 cm (9 1/2")



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Yasargil Micro Curette

Bayonet Shaped

Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Catalog #	Size
NS-635-01	24 cm (9 1/2")



Hardy Micro Sickle Dissector

Bayonet Shaped, Blunt


Catalog #	Size
NS-651-01	24 cm (9 1/2")
NS-651-02	24 cm (9 1/2")



Hardy Micro Sickle Dissector

Bayonet Shaped, Sharp

Catalog #	Size
NS-651-03	24 cm (9 1/2")
NS-651-04	24 cm (9 1/2")



Hardy Implant Fork

Bayonet Shaped


Catalog #	Size
NS-655-24	24 cm (9 1/2")



Hardy-Fahlbusch Micro Dissector

Bayonet Shaped Sharp

Catalog #	Size
NS-659-01	24 cm (9 1/2")
NS-659-02	24 cm (9 1/2")



INCHES
MILLIMETERS



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Casper Micro Dissector

Set of 6 Ref: NS-671-10 up to NS-675-10

Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Catalog #	Size
NS-670-00	N/A



Casper Micro Dissector

Straight

Catalog #	Size, Tip Size
NS-671-10	24 cm (9 1/2"), 1.0 mm
NS-671-45	24 cm (9 1/2"), 4.5 mm



NS-671-10




NS-671-45


Casper Micro Dissector

Bayonet Shaped, Curved Down


Catalog #	Size, Tip Size
NS-673-10	24 cm (9 1/2"), 1.0 mm
NS-673-20	24 cm (9 1/2"), 2.0 mm
NS-673-45	24 cm (9 1/2"), 4.5 mm



NS-673-10



NS-673-20




NS-673-45

Casper Micro Dissector

Bayonet Shaped, Curved Up


Catalog #	Size, Tip Size
NS-675-10	24 cm (9 1/2"), 1.0 mm




Micro Dissector

Straight

Catalog #	Size
NS-679-10	19 cm (7 1/2")
NS-679-22	21.5 cm (8 1/2")



NS-679-10



NS-679-22

NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.


Manufacturers of single use and reusable instruments



Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Ray Micro Curette

Angled 45°, Horizontal



Catalog

Size, Diameter

NS-681-01

19 cm (7 1/2"), Ø 4.0 mm

Ray Micro Curette

Angled 45°, Vertical



Catalog


Size, Diameter

NS-681-02

19 cm (7 1/2"), Ø 4.0 mm

Hardy Micro Curette

Bayonet Shaped



NS-685-04



NS-685-06



Catalog

Size, Diameter

NS-685-04


26 cm (10 1/4"), Ø 4.0 mm

NS-685-06

26 cm (10 1/4"), Ø 6.0 mm

Nicola Micro Curette

Bayonet Shaped, Left



Catalog

Size, Diameter

NS-689-01

21.5 cm (8 1/2"), Ø 6.5 mm

Nicola Micro Curette

Bayonet Shaped, Right



Catalog

Size, Diameter

NS-689-02

21.5 cm (8 1/2"), Ø 6.5 mm



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Nicola Micro Curette


Bayonet Shaped, Angled Up



Catalog #	Size, Diameter
NS-689-03	21.5 cm (8 1/2"), Ø 6.5 mm

Hosobuchi Fork Spatula

Bayonet Shaped, Rigid



Catalog #	Size
NS-695-26	26 cm (10 1/4")

Landolt Micro Curette


Bayonet Shaped, Malleable



Catalog #	Size
NS-697-26	26 cm (10 1/4")

Landolt Micro Curette


Bayonet Shaped, Blunt, Malleable



Catalog #	Size
NS-699-26	26 cm (10 1/4")

Landolt Micro Raspatory

Bayonet Shaped, Malleable



Catalog #	Size
NS-701-26	26 cm (10 1/4")

NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Landolt Micro Enucleator

Bayonet Shaped


Catalog #	Size
NS-703-01	26 cm (10 1/4")
NS-703-02	26 cm (10 1/4")



Landolt-Reulen Micro Dissector

Bayonet Shaped, Rigid

Catalog #	Size, Tip Size
NS-711-26	26 cm (10 1/4"), 2.0 mm



Landolt-Reulen Micro Hook

Bayonet Shaped, Blunt, Rigid


Catalog #	Size, Tip Size
NS-713-26	26 cm (10 1/4"), 1.7 mm



Landolt-Reulen Micro Curette

Bayonet Shaped, Rigid



Catalog #	Size, Diameter
NS-715-26	26 cm (10 1/4"), Ø 2.5 mm



Landolt-Reulen Micro Raspatory

Bayonet Shaped, Rigid

Catalog #	Size, Tip Size
NS-717-22	26 cm (10 1/4"), 2.2 mm
NS-717-32	26 cm (10 1/4"), 2.2 mm





CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments




Fahlbusch Micro Curette

Bayonet Shaped, Malleable

Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti


Catalog #	Size, Diameter
NS-719-01	26 cm (10 1/4"), Ø 2.5 mm



Fahlbusch Micro Scoop

Bayonet Shaped, Malleable

Catalog #	Size
NS-719-02	26 cm (10 1/4")




Cushing-Landolt Specula

Catalog #	Size, Blade Size
NS-751-01	N/A, 70 x 15 mm
NS-751-02	N/A, 90 x 15 mm
NS-751-03	N/A, 110 x 15 mm



Landolt Spcula

Catalog #	Size
NS-755-21	21.5 cm (8 1/2")



Mirror

Catalog #	Diameter
NS-790-06	Ø 6.0 mm
NS-790-08	Ø 8.0 mm
NS-790-10	Ø 10.0 mm
NS-790-12	Ø 12.0 mm



NEUROSURGERY



CENTRAL SURGICAL CO.(PVT)LTD.
Manufacturers of single use and reusable instruments



Micro Neuro Instruments
Mikro Neuro Instrumente
Micro Neuro Instrumentos
Micro Neuro Instruments
Micro Neuro Strumenti

Mirror Handle

Bayonet Shaped

Catalog #	Size
NS-792-00	15 cm (6")

